## AB-scope , Beauty Tester MicroScope

### The green fluorescence & White microscope


## Overview

#### Optimized for green fluorescence


*Using blue LED's for excitation with a 420-500nm emission filter designed to observe green fluorescence.*

**1.3** megapixel sensor

#### 1.3 Megapixels Sensor

*Experience and capture highly detailed images with 1.3 megapixel resolution*

#### High Magnification - 60x Fixed

*Up to 200x of glorious magnification depending of the distance of the Dino-Lite from your object.*

#### Measurement and Calibration Capabilities

*Ability to measure on pictures and live images. Also calibrate the microscope for assuring precision*

#### USB 2.0 Connection

*Digitize your microscopy experience; allowing you to easily record and share what you are looking at*

### A. Specifications

Model	DM5412VV-SP Dino-Lite Premier
Interface	USB 2.0
Product Resolution	1.3M pixels.
Magnification Rate	10x - 200x
Sensor	Color CMOS
Frame Rate	30fps
Save Formats	Image: BMP ,GIF ,PNG ,MNG ,TIF ,TGA ,PCX ,WBMP ,JP2 ,JPC ,JPG ,PGX ,RAS ,PNM
Microtouch	Touch sensitive trigger on the microscope for taking pictures
LED Lighting	7 Blue LED with 420-500nm emission filter and 1 White LED
Measurement Function	Yes
Calibration Function	Yes
Operating System Supported	Windows 7 ,Vista ,XP MAC OS 10.4 or later
Unit Weight	90(g)
Unit Dimension	10.5cm (H) x 3.2cm (D)
Package Dimensions	16cm (L) x 16cm (W) x 6cm (H)

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
### Illustration of skin Veins

**N**skin  
Normal skin

The ideal skin with sufficient moisture

*Tendency: Skin with an adequate amount of moisture and sebum*

---

This type of skin has sufficient moisture

Age	Laddie	Youth	Middle Age	Senium
Picture of Brow	20	30	40	50
Picture of eye wrinkle	20	30	40	50
Picture of arm	20	30	40	50
	Skin is Smooth soft and flexible.	Skin surface appeared slightly ageing furrow.	Skin surface appeared visible ageing furrow.	Skin surface appeared cloudiness and dim ageing furrow, but without flexible.


**Illustration of skin types**

The texture of skin is clean, small hairpore, surface smooth and not oily.  
**More moist but less sebaceous**

Neutral skin

Dry skin

Skin condition : dry, no polish, small pore and texture disordered.  
**Less moist and sebaceous**

**Ideal skin**  
The texture of skin is clean, moist fine, but not oily, soft and elastic.

The texture of skin is deep, hairpore big, and surface more oily.  
**More moist and sebaceous**

Oily skin

Dehydration but oily skin


Skin condition: coarse, Large Pore and texture not distinct.  
**Less moist but more sebaceous**

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


With out polarization

## Delicate Skin Rough skin

The skin does not have a sufficient amount of sebum. The protective sebaceous film is damaged so the skin is not protected properly. Moisturizing function is declined so moisture easily evaporates from the skin. As a result, the skin is susceptible to external irritation. When the skin is very external irritation. When the skin is very sensitive, it tends to become inflamed.

## Dark spots freckles (pigmentation)

Melanin is produced excessively because of the influence of UV rays. Melanin pigment remains within the skin because the skin's turnover is not smooth.

## Sebum plugs sebum

## Acne(Pigmentation)

Excessive sebum and old horny tissue harden and clog pores. They darken when oxidized. As this worsens, acne (pimple) develops.

There is too much sebum and the skin is sticky. Excessive sebum, dirt, and horny tissue have plugged the pore and acne has developed as a result. If you leave it, it will become inflamed.

## Skin tanned by UV-rays

## Dry skin/wrinkles caused by dryness

The skin is damaged by UV – rays. Valleys are not very clear and texture runs in one direction. The skin lacks moisture and very dehydrated. Advanced dry skin causes fine lines, sagging, and dark spots

Hills are flattened and valleys are unclear. Skin texture is rough and runs in one direction. The skin is dry, prone to roughness, and flat. It lacks both to moisture and sebum. The skin in this condition tends to become rough.

## Eye corners wrinkles

Dark spots and wrinkles become apparent on the eye corners as you age. Areas around the eye corners tend to have fine lines because they are prone to dryness and people tend to forget to care for them. Intensive moisturizing care is effective to treat them.

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